

Chromogranin A / CHGA (Neuroendocrine Marker) Antibody

Mouse Monoclonal Antibody [Clone LK2H10]

Catalog No	Format	Size
1113-MSM2-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
1113-MSM2-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
1113-MSM2-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug


Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	

Product Details


Clone	LK2H10
Gene Name	CHGA
Immunogen	Human pheochromocytoma
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1 / Kappa
Mol. Weight of Antigen	68-75kDa
Cellular Localization	Cytoplasmic vesicle, Neuronal dense core vesicle, Secreted, Secretory vesicle
Species Reactivity	Human, Monkey, Mouse, Pig, Rat
Positive Control	bowel, pancreas or pheochromocytoma., PC12 cells. Adrenal gland, Thyroid

**Optimal dilution for a specific application should be determined.*


Product Images for Chromogranin A / CHGA (Neuroendocrine Marker) Antibody




Formalin-fixed, paraffin-embedded human Small Cell Lung Carcinoma stained with Chromogranin A Monoclonal Antibody (LK2H10)




Formalin-fixed, paraffin-embedded Mouse Pancreas stained with Chromogranin A Monoclonal Antibody (LK2H10)




SDS-PAGE Analysis of Purified Chromogranin A Monoclonal Antibody (LK2H10). Confirmation of Purity and Integrity of Antibody.



Western blot analysis of human Pancreas tissue lysate using CHGA Mouse Monoclonal Antibody (LK2H10).



SDS-PAGE Analysis of Purified Chromogranin-A Mouse Monoclonal Antibody (LK2H10). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human Adrenal Gland stained with Chromogranin A Monoclonal Antibody (LK2H10)

Specificity & Comments

Chromogranin A is present in neuroendocrine cells throughout the body, including the neuroendocrine cells of the large and small intestine, adrenal medulla and pancreatic islets. It is an excellent marker for carcinoid tumors, pheochromocytomas, paragangliomas, and other neuroendocrine tumors. Co-expression of chromogranin A and neuron specific enolase (NSE) is common in neuroendocrine neoplasms. Reportedly, co-expression of certain keratins and chromogranin indicates neuroendocrine lineage. The presence of strong anti-chromogranin staining and absence of anti-keratin staining should raise the possibility of paraganglioma. The co-expression of chromogranin and NSE is typical of neuroendocrine neoplasms. Most pituitary adenomas and prolactinomas readily express chromogranin.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Cardiovascular, Immunology