

## DBC2 / RHOBTB2 (Tumor Suppressor) Antibody

Mouse Monoclonal Antibody [Clone DBC2/3364]

| Catalog No       | Format  | Size   |
|------------------|---|--------|
| 23221-MSM4-P0    | Purified Ab with BSA and Azide at 200ug/ml    | 20 ug  |
| 23221-MSM4-P1    | Purified Ab with BSA and Azide at 200ug/ml    | 100 ug |
| 23221-MSM4-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |


| Applications               | Tested Dillution | Note  |
|----------------------------|------------------|---|
| Immunohistochemistry (IHC) | 1-2ug/ml         | 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes |

### Product Details


|                        |  |
|------------------------|--|
| Clone                  | DBC2/3364  |
| Gene Name              | RHOBTB2  |
| Immunogen              | Human recombinant RHOBTB2 protein fragment (around aa 554-604) (exact sequence is proprietary) |
| Host                   | Mouse  |
| Clonality              | Monoclonal   |
| Isotype / Light Chain  | IgG2c / Kappa  |
| Mol. Weight of Antigen | 83kDa  |
| Species Reactivity     | Human  |
| Positive Control       | Human heart muscle or brain tissue.  |

\*Optimal dilution for a specific application should be determined.


### Product Images for DBC2 / RHOBTB2 (Tumor Suppressor) Antibody




Formalin-fixed, paraffin-embedded human thyroid stained with DBC2 Mouse Monoclonal Antibody (DBC2/3364).



Analysis of Protein Array containing more than 19,000 full-length human proteins using DBC2 Mouse Monoclonal Antibody (DBC2/3364). Z- and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody (MAb) (in combination with a fluorescently-tagged anti-IgG secondary antibody) produces when binding to a particular protein on the HuProt™ array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProt™ are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a MAb to its intended target. A MAb is considered to be specific to its intended target, if the MAb has an S-score of at least 2.5. For example, if a MAb binds to protein X with a Z-score of 43 and to protein Y with a Z-score of 14, then the S-score for the binding of that MAb to protein X is equal to 29.



SDS-PAGE Analysis of Purified Rho-related BTB domain-containing protein 2 Mouse Monoclonal Antibody (DBC2/3364). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human thyroid stained with DBC2 Mouse Monoclonal Antibody (DBC2/3364).

### Specificity & Comments

The Rho subfamily of Ras-related GTPases controls multiple aspects of cell function, including cytoskeletal rearrangement, nuclear signaling and cell growth. DBC-2 (deleted in breast cancer 2 gene protein), also known as RHOBTB2 (Rho-related BTB domain-containing protein 2), is a 727 amino acid member of the RhoBTB subfamily of Rho GTPases. Members of the RhoBTB subfamily are evolutionarily conserved and are characterized by a proline-rich region, a GTPase domain and two tandem BTB repeats. Expressed ubiquitously with highest levels in neural tissue, heart, brain and fetal lung, DBC-2 contains two BTB (POZ) domains through which it may bind to and regulate the function of target proteins, such as CUL-3. Additionally, DBC-2 is thought to function as a regulator of cell cycle and apoptosis events. Under normal conditions, DBC-2 is thought to exhibit tumor suppressor activity. Mutations in the gene encoding DBC-2 are associated with breast cancer, suggesting that mutated DBC-2 may play a role in carcinogenesis.

### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

### Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

### Research Areas

Signal Transduction