

## BOB.1 (B-Cell Marker) Antibody

Mouse Monoclonal Antibody [Clone BOB1/2424]


Catalog No	Format	Size
5450-MSM4-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
5450-MSM4-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
5450-MSM4-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Flow Cytometry (Flow)	1-2ug/million cells	
Immunofluorescence (IF)	1-3ug/ml	
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	


Product Details	
<b>Clone</b>	BOB1/2424
<b>Gene Name</b>	POU2AF1
<b>Immunogen</b>	Recombinant fragment (around aa 148-255) of human BOB1 (POU2AF1) protein (exact sequence is proprietary)
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype / Light Chain</b>	IgG2b / Kappa
<b>Mol. Weight of Antigen</b>	35kDa
<b>Cellular Localization</b>	Nucleus
<b>Species Reactivity</b>	Human
<b>Positive Control</b>	Raji, Ramos cells. Spleen or Tonsil.

\*Optimal dilution for a specific application should be determined.


### Product Images for BOB.1 (B-Cell Marker) Antibody




Western Blot Analysis of Raji cell lysate using BOB1 Monospecific Mouse Monoclonal Antibody (BOB1/2424).




SDS-PAGE Analysis Purified BOB1 Monospecific Mouse Monoclonal Antibody (BOB1/2424). Confirmation of Purity and Integrity of Antibody.




Flow Cytometric Analysis of paraformaldehyde-fixed Raji cells. BOB1 Monospecific Mouse Monoclonal Antibody (BOB1/2424); followed by goat anti-Mouse IgG-CF488 (blue) and Isotype Control (Red).



Immunofluorescent staining of paraformaldehyde-fixed Raji cells using BOB1 Monospecific Mouse Monoclonal Antibody (BOB1/2424) followed by goat anti-Mouse IgG conjugated to CF488 (green). Nuclei are labeled with RedDot (red).



Analysis of Protein Array containing >19,000 full-length human proteins using BOB1 Monospecific Mouse Monoclonal Antibody (BOB1/2424) Z- and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody (MAb) (in combination with a fluorescently-tagged anti-IgG secondary antibody) produces when binding to a particular protein on the HuProt™ array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProt™ are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a MAb to its intended target. A MAb is considered to specific to its intended target, if the MAb has an S-score of at least 2.5. For example, if a MAb binds to protein X with a Z-score of 43 and to protein Y with a Z-score of 14, then the S-score for the binding of that MAb to protein X is equal to 29.



Formalin-fixed, paraffin-embedded human Spleen stained with BOB1 Monospecific Mouse Monoclonal Antibody (BOB1/2424).

## Specificity & Comments

BOB.1 expression in a variety of established B-cell lines, representing different stages of B-cell development, has suggested a constitutive, B-cell-specific expression pattern. LP cells in nodular lymphocyte predominant Hodgkin lymphoma, because they are germinal center-derived, are consistently immuno-positive for BOB.1. Conversely, only some cases of classical Hodgkin lymphoma show BOB.1 immuno-reactivity within the Hodgkin and Reed-Sternberg cells. Expression of BOB.1 has been reported in follicular center cell lymphoma, diffuse large B-cell lymphoma and some cases of acute myeloid leukemia. B-CLL, marginal zone lymphoma, and mantle cell lymphoma may show weak to moderate immunoreactivity.

## Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

## Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

## Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

## Research Areas

B Cell Markers, Cancer, Cardiovascular