

Recombinant Rb1 (Tumor Suppressor Protein) Antibody

Rabbit Monoclonal Antibody [Clone RB1/2313R]

Catalog No	Format	Size
5925-RBM3-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
5925-RBM3-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
5925-RBM3-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug


Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	

Product Details


Clone	RB1/2313R
Gene Name	RB1
Immunogen	Recombinant human full-length Rb1 protein
Host	Rabbit
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	105kDa
Cellular Localization	Nucleus
Species Reactivity	Human, Mouse
Positive Control	A431 or BT-20 cells. Breast or Colon carcinoma., HeLa, Human Intestine, Human Kidney, Human Lung, Human Skin, Human Spleen

*Optimal dilution for a specific application should be determined.


Product Images for Recombinant Rb1 (Tumor Suppressor Protein) Antibody




Western Blot Analysis of Human Intestine, Mouse Intestine, Rat Intestine, Hamster Intestine, Guinea pig Intestine, Rabbit Intestine, Human Kidney, Mouse Kidney, Rat Kidney, Hamster Kidney, Guinea pig Kidney, Rabbit Kidney, Human Lung, Mouse Lung, Rat Lung, Hamster Lung, Guinea pig Lung and Rabbit Lung tissue lysates using Retinoblastoma Recombinant Rabbit Monoclonal Antibody (RB1/2313R).



Western Blot Analysis of Human Intestine, Human Kidney, Human Lung, Human Spleen and Human Skin tissue lysates using Retinoblastoma Recombinant Rabbit Monoclonal Antibody (RB1/2313R).



SDS-PAGE Analysis of Purified Retinoblastoma-associated protein Recombinant Rabbit Monoclonal Antibody (RB1/2313R). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human Colon Carcinoma stained with Retinoblastoma Recombinant Rabbit Monoclonal Antibody (RB1/2313R).

Specificity & Comments

Recognizes a 105kDa phosphoprotein, identified as retinoblastoma (Rb) gene product. Its epitope is localized between aa 703-772. It shows no cross reaction with p107 or p130. It specifically stains the nuclei of BT-20 cells and primary human foreskin fibroblast (HFF) cells. It does not stain the Rb-negative BT549 cells. It reacts with the hyperphosphorylated as well as the un (under) phosphorylated form of the Rb protein. Retinoblastoma gene product plays a key role in cell cycle control. It has been identified as a tumor suppressor gene whose loss of its function leads to tumor development. It is widely expressed in a variety of human tissues including breast, esophageal, squamous cell and cervical carcinoma.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Bladder Cancer, Breast Cancer, Cancer, Cardiovascular, Infectious Disease, Lung Cancer, Nuclear Marker, Transcription Factors