

## Bcl-X (Apoptosis Marker) Antibody

Mouse Monoclonal Antibody [Clone SPM337]

| Catalog No       | Format                                        | Size   |
|------------------|-----------------------------------------------|--------|
| 598-MSM2XX-P0    | Purified Ab with BSA and Azide at 200ug/ml    | 20 ug  |
| 598-MSM2XX-P1    | Purified Ab with BSA and Azide at 200ug/ml    | 100 ug |
| 598-MSM2XX-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |



| Applications               | Tested Dillution    | Note                                                                                                                                                                              |
|----------------------------|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Flow Cytometry (Flow)      | 1-2ug/million cells |                                                                                                                                                                                   |
| Immunofluorescence (IF)    | 1-3ug/ml            |                                                                                                                                                                                   |
| Immunohistochemistry (IHC) | 1-2ug/ml            | 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes |
| Western Blot (WB)          | 2-4ug/ml            |                                                                                                                                                                                   |

### Product Details

|                               |                                                                                                                                                                                                                                                    |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clone</b>                  | SPM337                                                                                                                                                                                                                                             |
| <b>Gene Name</b>              | BCL2L1                                                                                                                                                                                                                                             |
| <b>Immunogen</b>              | A synthetic peptide, aa 3-14 (Cys-QSNRELVVDLS) of human Bcl-X protein                                                                                                                                                                              |
| <b>Host</b>                   | Mouse                                                                                                                                                                                                                                              |
| <b>Clonality</b>              | Monoclonal                                                                                                                                                                                                                                         |
| <b>Isotype / Light Chain</b>  | IgG2a / Kappa                                                                                                                                                                                                                                      |
| <b>Mol. Weight of Antigen</b> | 27kDa                                                                                                                                                                                                                                              |
| <b>Cellular Localization</b>  | Centrosome, Cytoplasm, Cytoplasmic vesicle, Cytoskeleton, Cytosol, Microtubule organizing center, Mitochondrion inner membrane, Mitochondrion matrix, Mitochondrion outer membrane, Nucleus membrane, Secretory vesicle, Synaptic vesicle membrane |
| <b>Species Reactivity</b>     | Human, Mouse, Pig, Rat                                                                                                                                                                                                                             |
| <b>Positive Control</b>       | HL-60 or HeLa cells. Reed-Sternberg cells in Hodgkin's lymphoma., Jurkat, K562                                                                                                                                                                     |

\*Optimal dilution for a specific application should be determined.

### Product Images for Bcl-X (Apoptosis Marker) Antibody



SDS-PAGE Analysis of Purified bcl-x Mouse Monoclonal Antibody (SPM337). Confirmation of Purity and Integrity of Antibody.

Formalin-fixed, paraffin-embedded Hodgkin's lymphoma stained with Bcl-x Mouse Monoclonal Antibody (SPM337).

### Specificity & Comments

Recognizes a protein of 27kDa, identified as the Bcl-X protein. This MAb shows no cross-reaction with Bcl-2 or Bax protein. Bcl-X has two isoforms, Bcl-XL (long), a 241 amino acid protein which suppresses cell death. And Bcl-XS (short), a 178 amino acid protein lacking a 63 amino acid domain which functions as a dominant inhibitor of Bcl-2. This MAb reacts with both Bcl-XS and Bcl-XL proteins.

---

### Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

---

### Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

---

### Research Areas

Apoptosis, Autophagy, Cardiovascular, Colon Cancer, Cytokine Signaling, Immunology, Infectious Disease, Lung Cancer, Mitochondria Marker, Signal Transduction

---

### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

---