

Tal1 Antibody

Mouse Monoclonal Antibody [Clone TAL1/3146]

Catalog No	Format	Size
6886-MSM3-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
6886-MSM3-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
6886-MSM3-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug


Applications	Tested Dillution	Note
--------------	------------------	------

Product Details

Clone	TAL1/3146
Gene Name	TAL1
Immunogen	Recombinant human full-length TAL1 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1 / Kappa
Mol. Weight of Antigen	42kDa
Cellular Localization	Nucleus
Species Reactivity	Human
Positive Control	A549 xenograft.

**Optimal dilution for a specific application should be determined.*

Product Images for Tal1 Antibody



SDS-PAGE Analysis of Purified T-cell acute lymphocytic leukemia protein 1 Mouse Monoclonal Antibody (TAL1/3146). Confirmation of Purity and Integrity of Antibody.

Specificity & Comments

Activation of TAL1 characterizes up to 60% of cases of human T cell acute lymphoblastic leukemia, making it the most frequent gain-of-function mutation observed in this disorder. TAL1 (also designated SCL) is a serine phosphoprotein and basic helix-loop-helix transcription factor known to regulate embryonic hematopoiesis. This transcription factor binds as a heterodimer with E2A and HEB/HTF4 to a nucleotide sequence motif termed the E-box. In addition, leukemogenesis is accelerated dramatically by transgenic co-expression of TAL1 and the catalytic subunit of casein kinase II, a serine/threonine protein kinase known to modulate the activity of other HLH transcription factors.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Cardiovascular, Developmental Biology, Nuclear Marker, Transcription Factors