

CD71 / Transferrin Receptor (TFRC) Antibody

Mouse Monoclonal Antibody [Clone TFRC/1149]


| Catalog No | Format | Size |
|-----------------|-----------------------------------------------|--------|
| 7037-MSM4-P0 | Purified Ab with BSA and Azide at 200ug/ml | 20 ug |
| 7037-MSM4-P1 | Purified Ab with BSA and Azide at 200ug/ml | 100 ug |
| 7037-MSM4-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |

| Applications | Tested Dillution | Note |
|-------------------|------------------|------|
| Western Blot (WB) | 2-4ug/ml | |


| Product Details | |
|-------------------------------|------------------------------------------------|
| Clone | TFRC/1149 |
| Gene Name | TFRC |
| Immunogen | Recombinant human TFRC protein |
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype / Light Chain | IgG1 / Kappa |
| Mol. Weight of Antigen | 85-95kDa (monomer); 190kDa (dimer) |
| Cellular Localization | Cell membrane, Melanosome, Secreted |
| Species Reactivity | Human |
| Positive Control | PC3 cells. Human placenta or breast carcinoma. |

*Optimal dilution for a specific application should be determined.

Product Images for CD71 / Transferrin Receptor (TFRC) Antibody



SDS-PAGE Analysis of Purified Transferrin receptor protein 1 Mouse Monoclonal Antibody (TFRC/1149). Confirmation of Purity and Integrity of Antibody.



Western Blot of PC3 cell lysate using CD71 Monoclonal Antibody (TFRC/1149)

Specificity & Comments

It recognizes a ~90-95kDa protein which is identified as cell surface transferrin receptor (CD71), a disulfide-bonded homodimeric glycoprotein of 180-190kDa. This MAb is highly specific to CD71 and shows no cross-reaction with other related proteins. Ligand for transferrin receptor is the serum iron transport protein, transferrin. This receptor is broadly distributed in carcinomas, sarcomas, leukemias, and lymphomas. CD71/Transferrin receptor has been reported to be associated with cell proliferation in both normal and neoplastic tissues and useful in predicting clinical behavior or response to therapy in a number of malignancies including breast cancer.

Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Cancer, Cardiovascular, Signal Transduction

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.
