



Recombinant CD68 (Macrophage Marker) Antibody

Rabbit Monoclonal Antibody [Clone C68/3709R]

| Catalog No | Format | Size |
|-----------------|-----------------------------------------------|--------|
| 968-RBM12-P0 | Purified Ab with BSA and Azide at 200ug/ml | 20 ug |
| 968-RBM12-P1 | Purified Ab with BSA and Azide at 200ug/ml | 100 ug |
| 968-RBM12-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |

| Applications | Tested Dillution | Note |
|----------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunohistochemistry (IHC) | 1-2ug/ml | 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes |

Product Details

| | |
|------------------------|-------------------------------------------------------------------------------------|
| Clone | C68/3709R |
| Gene Name | CD68 |
| Immunogen | Recombinant fragment (around aa1-200) of human CD68 (exact sequence is proprietary) |
| Host | Rabbit |
| Clonality | Monoclonal |
| Isotype / Light Chain | IgG / Kappa |
| Mol. Weight of Antigen | 110kDa |
| Cellular Localization | Cell membrane, Endosome membrane, Lysosome membrane |
| Species Reactivity | Human |
| Positive Control | Human tonsil, lymph node or spleen. |

*Optimal dilution for a specific application should be determined.

Product Images for Recombinant CD68 (Macrophage Marker) Antibody

Specificity & Comments

This antibody recognizes a glycoprotein of 110kDa, which is identified as CD68. It is important for identifying macrophages in tissue sections. It stains macrophages in a wide variety of human tissues, including Kupffer cells and macrophages in the red pulp of the spleen, in lamina propria of the gut, in lung alveoli, and in bone marrow. It reacts with myeloid precursors and peripheral blood granulocytes. It also reacts with plasmacytoid T cells, which are supposed to be of monocyte/macrophage origin. It shows strong granular cytoplasmic staining of chronic and acute myeloid leukemia and also reacts with rare cases of true histiocytic neoplasia. Lymphomas are negative or show few granules.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Hematopoietic Stem Cells, Immunology