

CD79a (B-Cell Marker) Antibody

Mouse Monoclonal Antibody [Clone IGA/1406]

| Catalog No | Format | Size |
|----------------|---|--------|
| 973-MSM6-P0 | Purified Ab with BSA and Azide at 200ug/ml | 20 ug |
| 973-MSM6-P1 | Purified Ab with BSA and Azide at 200ug/ml | 100 ug |
| 973-MSM6-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |



| Applications | Tested Dillution | Note |
|-------------------------|---------------------|------|
| Flow Cytometry (Flow) | 1-2ug/million cells | |
| Immunofluorescence (IF) | 1-3ug/ml | |
| Western Blot (WB) | 2-4ug/ml | |

Product Details

| | |
|-------------------------------|---|
| Clone | IGA/1406 |
| Gene Name | CD79A |
| Immunogen | A synthetic peptide corresponding to aa 202-216 (GTYQDVGSLNIADVQ) of human CD79a protein. |
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype / Light Chain | IgG3 / Kappa |
| Mol. Weight of Antigen | 25-44kDa |
| Cellular Localization | Cell membrane |
| Species Reactivity | Human |
| Positive Control | Daudi or Ramos cells. Germinal center B- cells in a lymph node or tonsil. |

*Optimal dilution for a specific application should be determined.

Product Images for CD79a (B-Cell Marker) Antibody



SDS-PAGE Analysis of Purified CD79a Monoclonal Antibody (IGA/1406).
Confirmation of Integrity and Purity of Antibody.

Western Blot of Ramos cell lysate using CD79a Monoclonal Antibody (IGA/1406).

Specificity & Comments

A disulphide-linked heterodimer, consisting of mb-1 (or CD79a) and B29 (or CD79b) polypeptides, is non-covalently associated with membrane-bound immunoglobulins on B cells. This complex of mb-1 and B29 polypeptides and immunoglobulin constitute the B cell Ag receptor. CD79a first appears at pre B cell stage, early in maturation, and persists until the plasma cell stage where it is found as an intracellular component. CD79a is found in the majority of acute leukemias of precursor B cell type, in B cell lines, B cell lymphomas, and in some myelomas. It is not present in myeloid or T cell lines. Anti-CD79a is generally used to complement anti-CD20 especially for mature B-cell lymphomas after treatment with Rituximab (anti-CD20). This antibody will stain many of the same lymphomas as anti-CD20, but also is more likely to stain B-lymphoblastic lymphoma/leukemia than is anti-CD20. Anti-CD79a also stains more cases of plasma cell myeloma and occasionally some types of endothelial cells as well.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

B Cell Markers, Hematopoietic Stem Cells, Immunology, Infectious Disease, Mesenchymal Stem Cell Differentiation
