

Melanoma Marker (MART-1 + gp100) Antibody

Mouse Monoclonal Antibody [Clone DT101+ BC199 + HMB45]

Catalog No	Format	Size
MSM3-703-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
MSM3-703-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
MSM3-703-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug


Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details


Clone	DT101+ BC199 + HMB45
Gene Name	MLANA
Immunogen	Recombinant hMART-1 protein (DT101 & BC199); Extract of pigmented melanoma metastases from lymph nodes (HMB45)
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1, IgG2b / Kappa, Lambda
Mol. Weight of Antigen	20-22kDa (doublet) (MART-1): 90-100kDa (gp100)
Cellular Localization	Endoplasmic reticulum membrane, Golgi apparatus, Melanosome, trans-Golgi network membrane
Species Reactivity	Human
Positive Control	SK-MEL-13 and SK-MEL-19 melanoma cell lines. Human melanoma tissues.

*Optimal dilution for a specific application should be determined.

Product Images for Melanoma Marker (MART-1 + gp100) Antibody



Formalin-fixed, paraffin-embedded human Melanoma stained with Melanoma Marker Mouse Monoclonal Antibody (DT101+ BC199 + HMB45).



Formalin-fixed, paraffin-embedded human Melanoma stained with Melanoma Marker Mouse Monoclonal Antibody (DT101+ BC199 + HMB45).

Specificity & Comments

This antibody cocktail recognizes two melanoma-specific proteins, which include MART-1 and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Function of gp100 is not known but it is reported to be a useful marker for melanocytes and melanomas. A cocktail of these two markers is sensitive for labeling formalin-fixed, paraffin-embedded melanomas and other tumors showing melanocytic differentiation.

Supplied As

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

2 Union Square, Union City, CA 94587, Tel: 1-866-2-NEOBIO (1-866-263-6246)
orders@NeoBiotechnologies.com, www.NeoBiotechnologies.com