

BRCA1 (Breast Marker) Antibody

Mouse Monoclonal Antibody [Clone BRCA1/2973]

| Catalog No | Format | Size |
|----------------|---|--------|
| 672-MSM3-P0 | Purified Ab with BSA and Azide at 200ug/ml | 20 ug |
| 672-MSM3-P1 | Purified Ab with BSA and Azide at 200ug/ml | 100 ug |
| 672-MSM3-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |



| Applications | Tested Dillution | Note |
|--------------|------------------|------|
|--------------|------------------|------|

Product Details

| | |
|-------------------------------|---|
| Clone | BRCA1/2973 |
| Gene Name | BRCA1 |
| Immunogen | Recombinant fragment (around aa 445-620) of human BRCA1 protein (exact sequence is proprietary) |
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype / Light Chain | IgG2b / Lambda |
| Mol. Weight of Antigen | 220kDa |
| Cellular Localization | Cytoplasm. Nucleus. |
| Species Reactivity | Human |
| Positive Control | A431 or HeLa cells. Human breast ovarian carcinoma or renal cell carcinoma |

*Optimal dilution for a specific application should be determined.

Product Images for BRCA1 (Breast Marker) Antibody



SDS-PAGE Analysis of Purified BRCA1 Mouse Monoclonal Antibody (BRCA1/2973). Confirmation of Integrity and Purity of Antibody.

Analysis of Protein Array containing more than 19,000 full-length human proteins using BRCA-1 Mouse Monoclonal Antibody (BRCA1/2973) Z- and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody (MAb) (in combination with a fluorescently-tagged anti-IgG secondary antibody) produces when binding to a particular protein on the HuProt™ array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProt™ are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a MAb to its intended target. A MAb is considered to specific to its intended target, if the MAb has an S-score of at least 2.5. For example, if a MAb binds to protein X with a Z-score of 43 and to protein Y with a Z-score of 14, then the S-score for the binding of that MAb to protein X is equal to 29.

Specificity & Comments

This gene encodes a nuclear phosphoprotein that plays a role in maintaining genomic stability, and it also acts as a tumor suppressor. The encoded protein combines with other tumor suppressors, DNA damage sensors, and signal transducers to form a large multi-subunit protein complex known as the BRCA1-associated genome surveillance complex (BASC). This gene product associates with RNA polymerase II, and through the C-terminal domain, also interacts with histone deacetylase complexes. This protein thus plays a role in transcription, DNA repair of double-stranded breaks, and recombination. Mutations in this gene are responsible for approximately 40% of inherited breast cancers and more than 80% of inherited breast and ovarian cancers. Alternative splicing plays a role in modulating the subcellular localization and physiological function of this gene.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Breast Cancer, Cardiovascular, Infectious Disease, Nuclear Marker, Transcription Factors
