

## Recombinant Calbindin-1 / CALB1 (Neuronal Marker) Antibody

Mouse Monoclonal Antibody [Clone rKR6]

Catalog No	Format	Size
793-MSM10-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
793-MSM10-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
793-MSM10-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug


Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	

### Product Details


<b>Clone</b>	rKR6
<b>Gene Name</b>	CALB1
<b>Immunogen</b>	Prokaryotic recombinant protein corresponding to the majority of the calbindin molecule
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype / Light Chain</b>	IgG1 / Kappa
<b>Mol. Weight of Antigen</b>	28kDa
<b>Cellular Localization</b>	Cytoplasm
<b>Species Reactivity</b>	Human
<b>Positive Control</b>	Human brain, lung, pancreas or liver.

\*Optimal dilution for a specific application should be determined.


### Product Images for Recombinant Calbindin-1 / CALB1 (Neuronal Marker) Antibody



Western blot analysis of human brain tissue lysate using Calbindin 1 Recombinant Mouse Monoclonal Antibody (rKR6).



SDS-PAGE Analysis of Purified Calbindin Recombinant Mouse Monoclonal Antibody (rKR6). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human kidney stained with Calbindin Recombinant Mouse Monoclonal Antibody (rKR6). Inset: PBS instead of primary antibody; secondary only negative control.

---

### Specificity & Comments

Calbindin is a calcium-binding protein belonging to the troponin C superfamily. It functions as a buffer for cytoplasmic calcium and is found in the brain, kidney, gut and pancreatic islets. In normal brain, calbindin has been identified in medium sized neurons of the neuropil of the matrix compartment of the striatum, the woolly fiber arrangements of the globus pallidus and the fiber structures of the pars reticula of the substantia nigra. A proportion of dendritic cells within the light zone of germinal centers have also been noted to be positive for calbindin.

---

### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

---

### Supplied As

200ug/ml of Ab produced in CHO cell mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

---

### Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

---

### Research Areas

Neuroscience

---