

Recombinant IgD (Immunoglobulin Delta Heavy Chain) (B-Cell Marker) Antibody

Rabbit Monoclonal Antibody [Clone MSVA-701R]

| Catalog No | Format | Size |
|--------------|--------------------------------|--------|
| 3495-RBM9-P0 | Purified Ab with BSA and Azide | 20 ug |
| 3495-RBM9-P1 | Purified Ab with BSA and Azide | 100 ug |


| Applications | Tested Dillution | Note |
|----------------------------|------------------|--|
| Immunohistochemistry (IHC) | 1:100-1:200 | Manual Protocol: Freshly cut sections should be used (less than 10 days between cutting and staining). Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Target Retrieval Solution buffer. Apply the antibody at a dilution of 1:150 at 37°C for 60 minutes. Visualization of bound antibody by the EnVision Kit (Dako, Agilent) according to the manufacturer's directions. |

Product Details

| | |
|-------------------------------|---|
| Clone | MSVA-701R |
| Immunogen | Recombinant fragment (around aa1-200) of human IGHD protein (exact sequence is proprietary) |
| Host | Rabbit |
| Clonality | Monoclonal |
| Isotype / Light Chain | IgG / Kappa |
| Mol. Weight of Antigen | 44-80kDa (heavy chain); 21-25kDa (kappa light chain) |
| Cellular Localization | Cytoplasm |
| Species Reactivity | Human |
| Positive Control | Tonsil: A strong IgD staining should be seen in lymphocytes of the mantle zone. |

*Optimal dilution for a specific application should be determined.

Product Images for Recombinant IgD (Immunoglobulin Delta Heavy Chain) (B-Cell Marker) Antibody



Immunoglobulin heavy constant delta Rabbit Recombinant Monoclonal Antibody (MSVA-701R) tested on many normal and cancer tissues. The immunohistochemistry staining in these tissues aligns with the expression data in Human Protein Atlas.

Specificity & Comments

Immunoglobulins are four-chain, Y-shaped, monomeric structures comprised of two identical heavy chains and two identical light chains held together through interchain disulfide bonds. The chains form two domains, the Fab (antigen binding) fragment and the Fc (constant) fragment. Immunoglobulin D (IgD) exists as a monomer with 1 heavy chains and either 1 or 1 light chains. It plays a biological role as a transmembrane receptor molecule, co-expressed with IgM on the surface of mature/naive B cells. In particular, it is found on spleen B cell surfaces. Compared to IgM, IgD exists in much lower numbers and is not expressed on immature B cells. IgD surface expression on B cells is regulated in part by IL-27. In mice, the inhibition of this immunoglobulin isotype does not cause a significant change to the immune system.

Supplied As

Ab produced in HEK293 cell mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

B Cell Markers

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.
