

Recombinant Tyrosine Hydroxylase (Marker of Peripheral Neuroblastic Tumors) Antibody

Monoclonal Antibody [Clone TH/16182]

Catalog No	Format	Size
7054-MSM19-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
7054-MSM19-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
7054-MSM19-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug


Applications	Tested Dillution	Note
--------------	------------------	------

Product Details

Clone	TH/16182
Immunogen	Synthetic peptide corresponding to the C-terminal of human TH (epitope: aa520-528)
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	58kDa
Cellular Localization	Cell projection, Cytoplasm, Cytoplasmic vesicle, Nucleus
Species Reactivity	Human
Positive Control	kidney

**Optimal dilution for a specific application should be determined.*

Product Images for Recombinant Tyrosine Hydroxylase (Marker of Peripheral Neuroblastic Tumors) Antibody



SDS-PAGE Analysis of Purified Tyrosine 3-monooxygenase / TH Mouse Monoclonal Antibody (TH/16182). Confirmation of Purity and Integrity of Antibody.

Specificity & Comments

Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Neuroscience

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.