

## Recombinant GM-CSF (Granulocyte/Macrophage – Colony Stimulating Factor) Antibody

Rabbit Monoclonal Antibody [Clone CSF2/16944R]

Catalog No	Format	Size
1437-RBM11-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
1437-RBM11-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
1437-RBM11-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details	
Clone	CSF2/16944R
Immunogen	Recombinant full-length human CSF2 (Colony Stimulating Factor 2) protein
Host	Rabbit
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	16.29kDa
Cellular Localization	Secreted
Species Reactivity	Human

\*Optimal dilution for a specific application should be determined.

### Product Images for Recombinant GM-CSF (Granulocyte/Macrophage – Colony Stimulating Factor) Antibody

#### Specificity & Comments

Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.

#### Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.