

Recombinant CELA3B / ELA3B (Pancreatic Function Marker) Antibody

Mouse Monoclonal Antibody [Clone rCELA3B/1757]

Catalog No	Format	Size
23436-MSM15-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
23436-MSM15-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
23436-MSM15-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note

Product Details

Clone	rCELA3B/1757
Immunogen	Recombinant full-length human CELA3B protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1 / Kappa
Mol. Weight of Antigen	29.26kDa
Species Reactivity	Human
Positive Control	Pancreas

*Optimal dilution for a specific application should be determined.

Product Images for Recombinant CELA3B / ELA3B (Pancreatic Function Marker) Antibody

Specificity & Comments

Efficient protease with alanine specificity but only little elastolytic activity.

Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.