

Recombinant Mammaglobin (SCGB2A2) (Breast Cancer Marker) Antibody

Mouse Monoclonal Antibody [Clone rMGB/6618]

Catalog No	Format	Size
4250-MSM13-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
4250-MSM13-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
4250-MSM13-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
--------------	------------------	------

Product Details

Clone	rMGB/6618
Immunogen	Recombinant fragment of the human Mammaglobin protein (exact sequence is proprietary)
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1 / Kappa
Mol. Weight of Antigen	10.5kDa
Cellular Localization	Secreted
Species Reactivity	Human
Positive Control	Mammary gland specific

**Optimal dilution for a specific application should be determined.*

Product Images for Recombinant Mammaglobin (SCGB2A2) (Breast Cancer Marker) Antibody

Specificity & Comments

This antibody recognizes Mammaglobin-A. Clone: rMGB/6618.
Host: Mouse. Species reactivity: Human. Applications: ELISA.

Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.