

A1BG (Alpha-1-B glycoprotein) Antibody

Mouse Monoclonal Antibody [Clone A1BG/12985]

Catalog No	Format	Size
1-MSM9-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
1-MSM9-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
1-MSM9-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Western Blot (WB)	2-4ug/ml	

Product Details	
Clone	A1BG/12985
Immunogen	Recombinant human A1BG protein (exact sequence is proprietary)
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG
Mol. Weight of Antigen	54.25kDa
Cellular Localization	Secreted
Species Reactivity	Human
Positive Control	Plasma

*Optimal dilution for a specific application should be determined.

Product Images for A1BG (Alpha-1-B glycoprotein) Antibody

Specificity & Comments

This antibody recognizes Alpha-1B-glycoprotein. Clone: A1BG/12985. Host: Mouse. Applications: WB.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.