

Anti-mouse Cadherin-2 / CDH2 Antibody

Rat Monoclonal Antibody [Clone CDH2/13882]

Catalog No	Format	Size
1000-RTM24-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
1000-RTM24-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
1000-RTM24-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Western Blot (WB)	2-4ug/ml	

Product Details

Clone	CDH2/13882
Immunogen	Recombinant fragment (aa 308–597) corresponding to the extracellular domain of CDH2 protein
Host	Rat
Clonality	Monoclonal
Isotype / Light Chain	IgG2a / Kappa
Mol. Weight of Antigen	99.81kDa
Cellular Localization	Adherens junction, Cell junction, Cell membrane, Cell surface, Desmosome, Sarcolemma
Species Reactivity	Mouse

*Optimal dilution for a specific application should be determined.

Product Images for Anti-mouse Cadherin-2 / CDH2 Antibody

Specificity & Comments

Calcium-dependent cell adhesion protein; preferentially mediates homotypic cell-cell adhesion by dimerization with a CDH2 chain from another cell. Cadherins may thus contribute to the sorting of heterogeneous cell types. Acts as a regulator of neural stem cells quiescence by mediating anchorage of neural stem cells to ependymocytes in the adult subependymal zone: upon cleavage by MMP24, CDH2-mediated anchorage is affected, leading to modulate neural stem cell quiescence. Plays a role in cell-to-cell junction formation between pancreatic beta cells and neural crest stem (NCS) cells, promoting the formation of processes by NCS cells (By similarity). Required for proper neurite branching. Required for pre- and postsynaptic organization (By similarity). CDH2 may be involved in neuronal recognition mechanism. In hippocampal neurons, may regulate dendritic spine density.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.