

## Recombinant Mesothelin (Mesothelial Marker) Antibody

Rabbit Monoclonal Antibody [Clone MSLN/12899R]

Catalog No	Format	Size
10232-RBM21-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
10232-RBM21-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
10232-RBM21-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	

### Product Details

<b>Clone</b>	MSLN/12899R
<b>Immunogen</b>	Recombinant fragment (around aa 434-450) of human MSLN protein
<b>Host</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype / Light Chain</b>	IgG / Kappa
<b>Mol. Weight of Antigen</b>	68.98kDa
<b>Cellular Localization</b>	Cell membrane, Golgi apparatus, Secreted
<b>Species Reactivity</b>	Human
<b>Positive Control</b>	Expressed in lung

\*Optimal dilution for a specific application should be determined.

### Product Images for Recombinant Mesothelin (Mesothelial Marker) Antibody

#### Specificity & Comments

Membrane-anchored forms may play a role in cellular adhesion., Megakaryocyte-potentiating factor (MPF) potentiates megakaryocyte colony formation in vitro.

#### Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.