

## CSF2RA Antibody

Mouse Monoclonal Antibody [Clone CSF2RA/12584]

Catalog No	Format	Size
1438-MSM3-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
1438-MSM3-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
1438-MSM3-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Western Blot (WB)	2-4ug/ml	

### Product Details

Clone	CSF2RA/12584
Immunogen	Recombinant human CSF2RA protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2 / Kappa
Mol. Weight of Antigen	46.21kDa
Cellular Localization	Cell membrane, Secreted
Species Reactivity	Human

*\*Optimal dilution for a specific application should be determined.*

### Product Images for CSF2RA Antibody

#### Specificity & Comments

Low affinity receptor for granulocyte-macrophage colony-stimulating factor. Transduces a signal that results in the proliferation, differentiation, and functional activation of hematopoietic cells.

#### Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.