

Indoleamine 2,3-dioxygenase 2 Antibody

Mouse Monoclonal Antibody [Clone IDO2/4905]

Catalog No	Format	Size
169355-MSM5-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
169355-MSM5-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
169355-MSM5-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	IDO2/4905
Immunogen	Recombinant human IDO2 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2 / Kappa
Mol. Weight of Antigen	45.42kDa
Species Reactivity	Human
Positive Control	Detected in liver, small intestine, spleen, placenta, thymus, lung, brain, kidney, and colon (PubMed:17671174)

*Optimal dilution for a specific application should be determined.

Product Images for Indoleamine 2,3-dioxygenase 2 Antibody

Specificity & Comments

Catalyzes the first and rate limiting step of the catabolism of the essential amino acid tryptophan along the kynurenine pathway (PubMed:17671174). Involved in immune regulation. May not play a significant role in tryptophan-related tumoral resistance (PubMed:25691885).

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.