

## Recombinant Alpha-fetoprotein Antibody

Rabbit Monoclonal Antibody [Clone AFP/13720R]

Catalog No	Format	Size
174-RBM17-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
174-RBM17-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
174-RBM17-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

### Product Details

Clone	AFP/13720R
Immunogen	Recombinant full-length human AFP proteinAFP protein
Host	Rabbit
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	68.68kDa
Cellular Localization	Secreted
Species Reactivity	Human
Positive Control	Plasma

\*Optimal dilution for a specific application should be determined.

### Product Images for Recombinant Alpha-fetoprotein Antibody

#### Specificity & Comments

Binds copper, nickel, and fatty acids as well as, and bilirubin less well than, serum albumin. Only a small percentage (less than 2%) of the human AFP shows estrogen-binding properties.

#### Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.