

Friend leukemia integration 1 transcription factor (FLI1) Antibody

Mouse Monoclonal Antibody [Clone FLI1/6785]

Catalog No	Format	Size
2313-MSM6-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
2313-MSM6-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
2313-MSM6-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Western Blot (WB)	2-4ug/ml	

Product Details

Clone	FLI1/6785
Immunogen	Recombinant human FLI1 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1 / Kappa
Mol. Weight of Antigen	50.98kDa
Cellular Localization	Nucleus
Species Reactivity	Human

*Optimal dilution for a specific application should be determined.

Product Images for Friend leukemia integration 1 transcription factor (FLI1) Antibody

Specificity & Comments

Sequence-specific transcriptional activator (PubMed:24100448, PubMed:26316623, PubMed:28255014). Recognizes the DNA sequence 5'-C[CA]GGAAGT-3'.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.