

Zinc finger protein 324A Antibody

Mouse Monoclonal Antibody [Clone PCRP-ZNF324-1B7]

Catalog No	Format	Size
25799-MSM1-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
25799-MSM1-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
25799-MSM1-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Flow Cytometry (Flow)	1-2ug/million cells	
Immunofluorescence (IF)	1-3ug/ml	

Product Details	
Clone	PCRP-ZNF324-1B7
Immunogen	Recombinant human ZNF324 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2 / Kappa
Mol. Weight of Antigen	61.1kDa
Cellular Localization	Nucleus
Species Reactivity	Human
Positive Control	Expressed at high levels in the spleen, thymus, and PBMC, at low levels in the prostate, ovary, small intestine, colon (mucosal lining), placenta, lung, and pancreas, and very weakly expressed in the liver and kidney

**Optimal dilution for a specific application should be determined.*

Product Images for Zinc finger protein 324A Antibody

Specificity & Comments

May be involved in transcriptional regulation. May be involved in regulation of cell proliferation.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.