



Opticin Antibody

Mouse Monoclonal Antibody [Clone OPTC/13635]

Catalog No	Format	Size
26254-MSM1-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
26254-MSM1-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
26254-MSM1-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	OPTC/13635
Immunogen	Recombinant human OPTC protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG
Mol. Weight of Antigen	37.26kDa
Cellular Localization	Extracellular matrix, Extracellular space, Secreted
Species Reactivity	Human
Positive Control	Expressed in cartilage and synovial membranes (at protein level) (PubMed:18164633, PubMed:23845380)

**Optimal dilution for a specific application should be determined.*

Product Images for Opticin Antibody

Specificity & Comments

Inhibits angiogenesis in the vitreous humor of the eye, and therefore represses neovascularization (By similarity). Binds collagen fibrils (By similarity). May be involved in collagen fiber organization via regulation of other members of the small leucine-rich repeat proteoglycan superfamily (By similarity).

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.