

TOX high mobility group box family member 3 Antibody

Mouse Monoclonal Antibody [Clone PCRP-TOX3-1A6]

Catalog No	Format	Size
27324-MSM6-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
27324-MSM6-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
27324-MSM6-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Flow Cytometry (Flow)	1-2ug/million cells	
Immunofluorescence (IF)	1-3ug/ml	

Product Details

Clone	PCRP-TOX3-1A6
Immunogen	Recombinant human TOX3 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG1 / Kappa
Mol. Weight of Antigen	63.34kDa
Cellular Localization	Nucleus
Species Reactivity	Human
Positive Control	Expressed mainly in epithelial cells

**Optimal dilution for a specific application should be determined.*

Product Images for TOX high mobility group box family member 3 Antibody

Specificity & Comments

Transcriptional coactivator of the p300/CBP-mediated transcription complex. Activates transactivation through cAMP response element (CRE) sites. Protects against cell death by inducing antiapoptotic and repressing pro-apoptotic transcripts. Stimulates transcription from the estrogen-responsive or BCL-2 promoters. Required for depolarization-induced transcription activation of the C-FOS promoter in neurons. Associates with chromatin to the estrogen-responsive C3 promoter region.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.