

Recombinant Annexin A4 Antibody

Rabbit Monoclonal Antibody [Clone ANXA4/13470R]

Catalog No	Format	Size
307-RBM1-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
307-RBM1-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
307-RBM1-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Flow Cytometry (Flow)	1-2ug/million cells	
Immunofluorescence (IF)	1-3ug/ml	
Western Blot (WB)	2-4ug/ml	

Product Details

Clone	ANXA4/13470R
Immunogen	Recombinant fragment (around aa 2-139) of human ANXA4 protein
Host	Rabbit
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	35.88kDa
Cellular Localization	Zymogen granule membrane
Species Reactivity	Human

**Optimal dilution for a specific application should be determined.*

Product Images for Recombinant Annexin A4 Antibody

Specificity & Comments

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.