

## Recombinant Keratin, type II cytoskeletal 7 Antibody

Mouse Monoclonal Antibody [Clone rKRT7/760]

Catalog No	Format	Size
3855-MSM18-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
3855-MSM18-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
3855-MSM18-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	

### Product Details

<b>Clone</b>	rKRT7/760
<b>Immunogen</b>	Recombinant full-length human KRT7 protein
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype / Light Chain</b>	IgG1 / Kappa
<b>Mol. Weight of Antigen</b>	51.39kDa
<b>Cellular Localization</b>	Cytoplasm
<b>Species Reactivity</b>	Human
<b>Positive Control</b>	Expressed in cultured epidermal, bronchial and mesothelial cells but absent in colon, ectocervix and liver

*\*Optimal dilution for a specific application should be determined.*

### Product Images for Recombinant Keratin, type II cytoskeletal 7 Antibody

#### Specificity & Comments

Blocks interferon-dependent interphase and stimulates DNA synthesis in cells. Involved in the translational regulation of the human papillomavirus type 16 E7 mRNA (HPV16 E7).

#### Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.