

Brain-specific homeobox protein homolog Antibody

Mouse Monoclonal Antibody [Clone PCR-P-BSX-5H1]

Catalog No	Format	Size
390259-MSM1-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
390259-MSM1-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
390259-MSM1-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Flow Cytometry (Flow)	1-2ug/million cells	
Immunofluorescence (IF)	1-3ug/ml	

Product Details

Clone	PCR-P-BSX-5H1
Immunogen	Recombinant human SX3 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG
Mol. Weight of Antigen	25.93kDa
Cellular Localization	Nucleus
Species Reactivity	Human

*Optimal dilution for a specific application should be determined.

Product Images for Brain-specific homeobox protein homolog Antibody

Specificity & Comments

DNA binding protein that function as transcriptional activator. Is essential for normal postnatal growth and nursing. Is an essential factor for neuronal neuropeptide Y and agouti-related peptide function and locomotory behavior in the control of energy balance (By similarity).

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.