

## Neurofilament (NF-H) (Neuronal Marker) Antibody

Mouse Monoclonal Antibody [Clone NEFH/15686]

Catalog No	Format	Size
4744-MSM8-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
4744-MSM8-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
4744-MSM8-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Western Blot (WB)	2-4ug/ml	

### Product Details

<b>Clone</b>	NEFH/15686
<b>Immunogen</b>	Recombinant fragment (around aa 327-468) of human NEFH protein
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype / Light Chain</b>	IgG1
<b>Mol. Weight of Antigen</b>	111.84kDa
<b>Cellular Localization</b>	Axon, Cell projection, Cytoplasm, Cytoskeleton
<b>Species Reactivity</b>	Human

*\*Optimal dilution for a specific application should be determined.*

### Product Images for Neurofilament (NF-H) (Neuronal Marker) Antibody

#### Specificity & Comments

Neurofilaments usually contain three intermediate filament proteins: NEFL, NEFM, and NEFH which are involved in the maintenance of neuronal caliber. NEFH has an important function in mature axons that is not subserved by the two smaller NF proteins. May additionally cooperate with the neuronal intermediate filament proteins PRPH and INA to form neuronal filamentous networks (By similarity).

#### Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.