

## NTRK2 (Neurotrophic Receptor Tyrosine Kinase 2) Antibody

Mouse Monoclonal Antibody [Clone NTRK2/6718]

| Catalog No      | Format                                       | Size   |
|-----------------|----------------------------------------------|--------|
| 4915-MSM1-P0    | Purified Ab with BSA and Azide at 200ug/ml   | 20 ug  |
| 4915-MSM1-P1    | Purified Ab with BSA and Azide at 200ug/ml   | 100 ug |
| 4915-MSM1-P1ABX | Purified Ab WITHOUT BSA or Azide at 1.0mg/ml | 100 ug |

| Applications               | Tested Dillution | Note                                                                                                                                                                              |
|----------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunohistochemistry (IHC) | 1-2ug/ml         | 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes |

### Product Details

|                               |                                                                                                                                                 |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clone</b>                  | NTRK2/6718                                                                                                                                      |
| <b>Immunogen</b>              | Recombinant human NTRK1 protein                                                                                                                 |
| <b>Host</b>                   | Mouse                                                                                                                                           |
| <b>Clonality</b>              | Monoclonal                                                                                                                                      |
| <b>Isotype / Light Chain</b>  | IgG1 / Kappa                                                                                                                                    |
| <b>Mol. Weight of Antigen</b> | 92.0kDa                                                                                                                                         |
| <b>Cellular Localization</b>  | Axon, Cell membrane, Cell projection, Cytoplasm, Dendrite, Early endosome membrane, Endosome membrane, Perinuclear region, Postsynaptic density |
| <b>Species Reactivity</b>     | Human                                                                                                                                           |
| <b>Positive Control</b>       | Isoform TrkB is expressed in the central and peripheral nervous system                                                                          |

\*Optimal dilution for a specific application should be determined.

### Product Images for NTRK2 (Neurotrophic Receptor Tyrosine Kinase 2) Antibody

#### Specificity & Comments

Receptor tyrosine kinase involved in the development and the maturation of the central and the peripheral nervous systems through regulation of neuron survival, proliferation, migration, differentiation, and synapse formation and plasticity (By similarity). Receptor for BDNF/brain-derived neurotrophic factor and NTF4/neurotrophin-4. Alternatively can also bind NTF3/neurotrophin-3 which is less efficient in activating the receptor but regulates neuron survival through NTRK2 (PubMed:15494731, PubMed:7574684). Upon ligand-binding, undergoes homodimerization, autophosphorylation and activation (PubMed:15494731). Recruits, phosphorylates and/or activates several downstream effectors including SHC1, FRS2, SH2B1, SH2B2 and PLCG1 that regulate distinct overlapping signaling cascades. Through SHC1, FRS2, SH2B1, SH2B2 activates the GRB2-Ras-MAPK cascade that regulates for instance neuronal differentiation including neurite outgrowth. Through the same effectors controls the Ras-PI3 kinase-AKT1 signaling cascade that mainly regulates growth and survival. Through PLCG1 and the downstream protein kinase C-regulated pathways controls synaptic plasticity. Thereby, plays a role in learning and memory by regulating both short term synaptic function and long-term potentiation. PLCG1 also leads to NF-Kappa-B activation and the transcription of genes involved in cell survival. Hence, it is able to suppress anoikis, the apoptosis resulting from loss of cell-matrix interactions. May also play a role in neutrophin-dependent calcium signaling in glial cells and mediate communication between neurons and glia.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

#### Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.