

Serpin B5 Antibody

Mouse Monoclonal Antibody [Clone SERPINB5/4976]

Catalog No	Format	Size
5268-MSM6-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
5268-MSM6-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
5268-MSM6-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	SERPINB5/4976
Immunogen	Recombinant fragment of human SERPINB5 protein (around aa1-200)
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2b / Kappa
Mol. Weight of Antigen	42.1kDa
Cellular Localization	Extracellular space, Secreted
Species Reactivity	Human
Positive Control	Normal mammary epithelial cells

*Optimal dilution for a specific application should be determined.

Product Images for Serpin B5 Antibody

Specificity & Comments

Tumor suppressor. It blocks the growth, invasion, and metastatic properties of mammary tumors. As it does not undergo the S (stressed) to R (relaxed) conformational transition characteristic of active serpins, it exhibits no serine protease inhibitory activity.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.