

Uroplakin 1B (Urothelial Differentiation Marker) Antibody

Mouse Monoclonal Antibody [Clone UPK1B/12693]

Catalog No	Format	Size
7348-MSM13-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
7348-MSM13-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
7348-MSM13-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	UPK1B/12693
Immunogen	Recombinant fragment (around aa 109-229) of human UPK1B protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2b / Kappa
Mol. Weight of Antigen	29.64kDa
Cellular Localization	Membrane
Species Reactivity	Human
Positive Control	Bladder epithelium

**Optimal dilution for a specific application should be determined.*

Product Images for Uroplakin 1B (Urothelial Differentiation Marker) Antibody

Specificity & Comments

Component of the asymmetric unit membrane (AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in normal bladder epithelial physiology, possibly in regulating membrane permeability of superficial umbrella cells or in stabilizing the apical membrane through AUM/cytoskeletal interactions (By similarity).

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.