

Regulator of G-protein signaling 5 Antibody

Mouse Monoclonal Antibody [Clone BICCN-RGS5-1A11]

Catalog No	Format	Size
8490-MSM3-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
8490-MSM3-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
8490-MSM3-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	BICCN-RGS5-1A11
Immunogen	Recombinant fragment (around aa 41-174) of human RGS5 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2a / Kappa
Mol. Weight of Antigen	20.95kDa
Cellular Localization	Cytoplasm, Membrane
Species Reactivity	Human

**Optimal dilution for a specific application should be determined.*

Product Images for Regulator of G-protein signaling 5 Antibody

Specificity & Comments

Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form. Binds to G(i)-alpha and G(o)-alpha, but not to G(s)-alpha (By similarity).

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.