



## Recombinant Cyclin E (G1/S-Phase Cyclin) Antibody

Mouse Monoclonal Antibody [Clone rCCNE1/13816]

Catalog No	Format	Size
898-MSM10-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
898-MSM10-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
898-MSM10-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes
Western Blot (WB)	2-4ug/ml	

### Product Details

Clone	rCCNE1/13816
Immunogen	Recombinant fragment (around aa 10-176) of human CCNE1 protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2b / Kappa
Mol. Weight of Antigen	47.08kDa
Cellular Localization	Nucleus
Species Reactivity	Human
Positive Control	Highly expressed in testis and placenta

*\*Optimal dilution for a specific application should be determined.*

### Product Images for Recombinant Cyclin E (G1/S-Phase Cyclin) Antibody

#### Specificity & Comments

Essential for the control of the cell cycle at the G1/S (start) transition.

#### Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

#### Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

#### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.