

Recombinant Interleukin-33 Antibody

Rabbit Monoclonal Antibody [Clone IL33/17275R]

Catalog No	Format	Size
90865-RBM2-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
90865-RBM2-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
90865-RBM2-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes


Product Details

Clone	IL33/17275R
Immunogen	Recombinant human IL33 protein
Host	Rabbit
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	30.76kDa
Cellular Localization	Chromosome, Cytoplasm, Cytoplasmic vesicle, Nucleus, Secreted, Secretory vesicle
Species Reactivity	Human
Positive Control	Expressed at high level in high endothelial venules found in tonsils, Peyer patches and mesenteric lymph nodes

*Optimal dilution for a specific application should be determined.

Product Images for Recombinant Interleukin-33 Antibody

Purity: SDS-PAGE



SDS-PAGE Analysis of Purified Recombinant Interleukin-33 Rabbit Monoclonal Antibody (IL33/17275R). Confirmation of Purity and Integrity of Antibody.

Specificity & Comments

Cytokine that binds to and signals through the IL1RL1/ST2 receptor which in turn activates NF-kappa-B and MAPK signaling pathways in target cells (PubMed:16286016, PubMed:19841166). Involved in the maturation of Th2 cells inducing the secretion of T-helper type 2-associated cytokines (PubMed:17853410, PubMed:18836528). Also involved in activation of mast cells, basophils, eosinophils and natural killer cells (PubMed:17853410, PubMed:18836528). Acts as an enhancer of polarization of alternatively activated macrophages (PubMed:19841166). Acts as a chemoattractant for Th2 cells, and may function as an 'alarmin', that amplifies immune responses during tissue injury (PubMed:17853410, PubMed:18836528). Induces rapid UCP2-dependent mitochondrial rewiring that attenuates the generation of reactive oxygen species and preserves the integrity of Krebs cycle required for persistent production of itaconate and subsequent GATA3-dependent differentiation of inflammation-resolving alternatively activated macrophages (By similarity)., In quiescent endothelia the uncleaved form is constitutively and abundantly expressed, and acts as a chromatin-associated nuclear factor with transcriptional repressor properties, it may sequester nuclear NF-kappaB/RELA, lowering expression of its targets (PubMed:21734074). This form is rapidly lost upon angiogenic or pro-inflammatory activation (PubMed:18787100).

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.