

L-lactate dehydrogenase A-like 6B Antibody

Mouse Monoclonal Antibody [Clone LDHAL6B/12744]

Catalog No	Format	Size
92483-MSM4-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
92483-MSM4-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
92483-MSM4-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Western Blot (WB)	2-4ug/ml	

Product Details

Clone	LDHAL6B/12744
Immunogen	Recombinant human LDHAL6B protein
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2b / Kappa
Mol. Weight of Antigen	41.94kDa
Species Reactivity	Human
Positive Control	Testis specific

**Optimal dilution for a specific application should be determined.*

Product Images for L-lactate dehydrogenase A-like 6B Antibody

Specificity & Comments

This antibody recognizes L-lactate dehydrogenase A-like 6B. Clone: LDHAL6B/12744. Host: Mouse. Applications: WB.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.