

C-X-C motif chemokine 14 Antibody

Mouse Monoclonal Antibody [Clone CXCL14/]

Catalog No	Format	Size
9547-MSM1-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
9547-MSM1-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
9547-MSM1-P1ABX	Purified Ab WITHOUT BSA or Azide at 1.0mg/ml	100 ug

Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	CXCL14/
Immunogen	Recombinant fragment (around aa 35-111) of human CXCL14 protein (exact sequence is proprietary)
Host	Mouse
Clonality	Monoclonal
Isotype / Light Chain	IgG2c / Kappa
Mol. Weight of Antigen	13.08kDa
Cellular Localization	Secreted
Species Reactivity	Human
Positive Control	Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas

*Optimal dilution for a specific application should be determined.

Product Images for C-X-C motif chemokine 14 Antibody

Specificity & Comments

Potent chemoattractant for neutrophils, and weaker for dendritic cells. Not chemotactic for T-cells, B-cells, monocytes, natural killer cells or granulocytes. Does not inhibit proliferation of myeloid progenitors in colony formation assays.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.