

Recombinant Pan-Cadherin Antibody

Rabbit Monoclonal Antibody [Clone Pan-CAD/12808R]

Catalog No	Format	Size
RBM4-12808-P0	Purified Ab with BSA and Azide at 200ug/ml	20 ug
RBM4-12808-P1	Purified Ab with BSA and Azide at 200ug/ml	100 ug
RBM4-12808-P1ABX	Purified Ab WITHOUT BSA and Azide at 1.0mg/ml	100 ug



Applications	Tested Dillution	Note
Immunohistochemistry (IHC)	1-2ug/ml	30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes

Product Details

Clone	Pan-CAD/12808R
Gene Name	CDH1
Immunogen	Recombinant fragment corresponding to the highly conserved C-terminal cytoplasmic tail consensus sequence (aa 892–906) of human classical cadherins
Host	Rabbit
Clonality	Monoclonal
Isotype / Light Chain	IgG / Kappa
Mol. Weight of Antigen	120 kDa
Cellular Localization	Cell membrane
Species Reactivity	Human
Positive Control	Human heart, colon, prostate or brain.

*Optimal dilution for a specific application should be determined.

Product Images for Recombinant Pan-Cadherin Antibody



SDS-PAGE Analysis of Purified pan-Cadherin Recombinant Rabbit Monoclonal Antibody (Pan-CAD/12808R). Confirmation of Integrity and Purity of Antibody.

Formalin-fixed, paraffin-embedded human colon stained with Pan Cadherin Recombinant Rabbit Monoclonal Antibody (Pan-CAD/12808R). HIER: Tris/EDTA, pH9.0, 45min. 2°: HRP-polymer, 30min. DAB, 5min.

Specificity & Comments

Cadherins comprise a family of Ca²⁺-dependent adhesion molecules that function to mediate cell-cell binding critical to the maintenance of tissue structure and morphogenesis. The classical cadherins, E-, N- and P-cadherin, consist of large extracellular domains characterized by a series of five homologous NH₂ terminal repeats. The most distal of these cadherins is thought to be responsible for binding specificity, transmembrane domains and carboxy terminal intracellular domains. The relatively short intracellular domains interact with a variety of cytoplasmic proteins, such as β -catenin, to regulate cadherin function. Members of this family of adhesion proteins include rat cadherin K (and its human homolog, cadherin 6), R-cadherin, B-cadherin, E/P cadherin and cadherin-5.

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

Supplied As

200ug/ml of Ab produced in a mammalian-based expression system. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Bladder Cancer, Colon Cancer, Immunology, Infectious Disease, Signal Transduction
